




At the heart of
critical moment


Central Venous Catheters



www.NSDmedical.com

 Office A1 - 327B, P.O.Box : 932
Ajman FZ Building, U.A.E

 +971 50 641 98 73

 www.nsdmedical.com
Nelson@nsdmedical.com

Invasive Blood pressure Transducer

- Integrated unit for clinical accuracy
- Dual Mode of flushing: Controlled squeeze & Snap-tab
- Completely transparent for ease of priming and debubbling
- Available with flushing of 3cc/hr, 30cc/hr or without flush
- Easy Zeroing



Ancillary Products

- Pressure Monitoring Lines
- Pediatric Pressure Monitoring Lines
- Low Prime Lines
- Extension Lines (with or without stopcock)
- Flow Regulator lines
- IV Filter Lines
- Flo Regulator
- MV Set
- Multiple Extension with Needleless Valves



Class VI implant tested tubing, DEHP Free
 Available in wide variants of Lengths, Internal Diameter, Wall Thickness and Colours



Contact for your exact requirement.


Flow Regulator	B-FR-01P
----------------	----------



Measure Volume Set	B-BR-150Y
--------------------	-----------



Needleless Valve with Multiple Extension	Needleless Valve Single Double Triple	BEXNV BEXNV-S BEXNV-D BEXNV-T
--	--	--




Pressure Monitoring Kits

- Choice of standard Single, Double and Triple kits
- Custom Design kit as per user specifications and requirement
- Option of Closed loop blood sampling kit available including needleless sampling site
- Class VI implant tested tubing, DEHP Free
- Contact for non-obligatory schematic design / samples



Schematic design of BKT-165 (Standard Double Kit)



Pressure Bag / Infuser

- Disposable 500ml B-PI-500D
- Disposable 1000ml B-PI-1000D
- Reusable 500ml B-PI-500R
- Reusable 1000ml B-PI-1000R



Cables available for all types of monitors

Ordering Information

	Catheter Size Ga/Fr	Lumen Ga	Cather Length cms	Dilator Fr	Guide Wire Nickle Titanium Double Distal	with Std. Needle	with Y Needle Std. Syringe	with Std. Needle Introducer Syringe
						P/N	P/N	P/N
Single Lumen Catheter								
1 Lumen 22Ga 4cms	22Ga	22	4	4	0.018x30 cm	NSD-CVC-1001	---	---
1 Lumen 22Ga 4cms with 22G Cannula	22Ga	22	4	4	0.018x45 cm	NSD-CVC-1002	---	---
1 Lumen 22Ga 8cms	22GA	22	8	4	0.018x45 cm	NSD-CVC-1003	---	---
1 Lumen 22Ga 10cms	22Ga	22	10	4	0.018x45 cm	NSD-CVC-1004	---	---
1 Lumen 22Ga 13cms	22 Ga	22	13	4	0.018x45 cm	NSD-CVC-1005	---	---
1 Lumen 22Ga 40cms	22 GA	22	40	4	0.018x45 cm	NSD-CVC-1006	---	---
1 Lumen 20Ga 13cms	20 Ga	20	13	4	0.018x45 cm	NSD-CVC-1007	---	---
1 Lumen 18Ga 10cms	18 Ga	18	10	5	0.018x45 cm	NSD-CVC-1008	NSD-CVC-18G01-10YN	NSD-CVC-18G01-10IS
1 Lumen 18Ga 13cms	18GA	18	13	5	0.018x45 cm	NSD-CVC-1009	NSD-CVC-18G01-13YN	NSD-CVC-18G01-13IS
1 Lumen 18Ga 20cms	18 Ga	18	20	5	0.018x45 cm	NSD-CVC-1010	NSD-CVC-18G01-20YN	NSD-CVC-18G01-20IS
1 Lumen 16Ga 13cms	16Ga	16	13	5	0.035x45 cm	NSD-CVC-1011	NSD-CVC-16G01-13YN	NSD-CVC-16G01-13IS
1 Lumen 16Ga 16cms	16 Ga	16	16	5	0.035x45 cm	NSD-CVC-1012	NSD-CVC-16G01-16YN	NSD-CVC-16G01-16IS
1 Lumen 16Ga 20cms	16 Ga	16	20	5	0.035x60 cm	NSD-CVC-1013	NSD-CVC-16G01-20YN	NSD-CVC-16G01-20IS
1 Lumen 16Ga 40cms	16GA	16	40	5	0.035x60 cm	NSD-CVC-1014	NSD-CVC-16G01-40YN	NSD-CVC-16G01-40IS
1 Lumen 14Ga 20cms	14 Ga	14	20	8	0.035x60 cm	NSD-CVC-1015	NSD-CVC-14G01-20YN	NSD-CVC-14G01-20IS
Double Lumen Catheter								
2 Lumen 4Fr 6cms	4 Fr	20, 22	6	5	0.018x 45 cm	NSD-CVC-2001	NSD-CVC-402-06YN	NSD-CVC-402-6IS
2 Lumen 4Fr 8cms	4 Fr	20, 22	8	5	0.018x 45 cm	NSD-CVC-2002	NSD-CVC-402-08YN	NSD-CVC-402-8IS
2 Lumen 4Fr 13cms	4 Fr	20, 22	13	5	0.018x 45 cm	NSD-CVC-2003	NSD-CVC-402-13YN	NSD-CVC-402-13IS
2 Lumen 5Fr 8cms	5 Fr	18, 20	8	6	0.018x 45 cm	NSD-CVC-2004	NSD-CVC-502-08YN	NSD-CVC-502-8IS
2 Lumen 5Fr 13cms	5 Fr	18, 20	13	6	0.018x 45 cm	NSD-CVC-2005	NSD-CVC-502-13YN	NSD-CVC-502-13IS
2 Lumen 7Fr 13cms	7 Fr	14, 18	13	8	0.035x 45 cm	NSD-CVC-2006	NSD-CVC-702-13YN	NSD-CVC-702-13IS
2 Lumen 7Fr 16cms	7 Fr	14, 18	16	8	0.035x 45 cm	NSD-CVC-2007	NSD-CVC-702-16YN	NSD-CVC-702-16IS
2 Lumen 7Fr 20cms	7 Fr	14, 18	20	8	0.035x 60 cm	NSD-CVC-2008	NSD-CVC-702-20YN	NSD-CVC-702-20IS
Triple Lumen Catheter								
3 Lumen 4.5Fr 6cms with 22G Cannula	4.5 Fr	20, 22, 22	6	5	0.018x45 cm	NSD-CVC-3001	---	---
3 Lumen 4.5Fr 6cms	4.5 Fr	20, 22, 22	6	5	0.018x 45 cm	NSD-CVC-3002	NSD-CVC-4.503-06YN	NSD-CVC-4.503-6IS
3 Lumen 4.5Fr 10cms	4.5 Fr	20, 22, 22	10	5	0.018x 45 cm	NSD-CVC-3003	NSD-CVC-4.503-10YN	NSD-CVC-4.503-10IS
3 Lumen 5.5Fr 6cms	5.5 Fr	20, 22, 22	6	6	0.018x 45 cm	NSD-CVC-3004	NSD-CVC-5.503-06YN	NSD-CVC-5.503-6IS
3 Lumen 5.5Fr 8cms	5.5 Fr	20, 22, 22	8	6	0.018x 45 cm	NSD-CVC-3005	NSD-CVC-5.503-08YN	NSD-CVC-5.503-8IS
3 Lumen 5.5Fr 13cms	5.5 Fr	20, 22, 22	13	6	0.018x 45 cm	NSD-CVC-3006	NSD-CVC-5.503-13YN	NSD-CVC-5.503-13IS
3 Lumen 7Fr 10cms	7 Fr	16, 18, 18	10	8	0.035x 45 cm	NSD-CVC-3007	NSD-CVC-703-10YN	NSD-CVC-703-10IS
3 Lumen 7Fr 13cms	7 Fr	16, 18, 18	13	8	0.035x 45 cm	NSD-CVC-3008	NSD-CVC-703-13YN	NSD-CVC-703-13IS
3 Lumen 7Fr 16cms	7 Fr	16, 18, 18	16	8	0.035x 45 cm	NSD-CVC-3009	NSD-CVC-703-16YN	NSD-CVC-703-16IS
3 Lumen 7Fr 20cms	7 Fr	16, 18, 18	20	8	0.035x 60 cm	NSD-CVC-3010	NSD-CVC-703-20YN	NSD-CVC-703-20IS
Four Lumen Catheter								
4 Lumen 8.5Fr 10cms	8.5 Fr	14,16,18,18	10	9	0.035x 45 cm	NSD-CVC-4001	NSD-CVC-8.504-10YN	NSD-CVC-8.504-10IS
4 Lumen 8.5Fr 13cms	8.5 Fr	14,16,18,18	13	9	0.035x 45 cm	NSD-CVC-4002	NSD-CVC-8.504-13YN	NSD-CVC-8.504-13IS
4 Lumen 8.5Fr 16cms	8.5 Fr	14,16,18,18	16	9	0.035x 45 cm	NSD-CVC-4003	NSD-CVC-8.504-16YN	NSD-CVC-8.504-16IS
4 Lumen 8.5Fr 20cms	8.5 Fr	14,16,18,18	20	9	0.035x 60 cm	NSD-CVC-4004	NSD-CVC-8.504-20YN	NSD-CVC-8.504-20IS

Option of Needle-less valve is also available in place of Injection port



18 St James Road, Purley, Surrey,
Greater London, CR 8 2DL, United Kingdom
+44 207 305 7173

www.nsdmedical.com
Sales@nsdmedical.com

NSD/BL/CAT/04.
Rev05 19/9/2023